

# O<sub>2</sub> Fingerprinting in Pasteurised Beverage Products Using Headspace Gas Analysis (Orbisphere 6110)

## Objective

This report aims to evaluate the use of the Orbisphere 6110 for headspace gas analysis in pasteurised carbonated beverage packaging. The specific goal is to develop an O<sub>2</sub> fingerprinting approach to identify post-pasteurisation oxygen ingress. A particular challenge is that dissolved oxygen can react with the beer matrix over time, rendering direct O<sub>2</sub> measurement increasingly unreliable. Therefore, a derived metric—**residual gas volume** adjusted for CO<sub>2</sub> purity—can be used to infer the **historical presence and consumption** of oxygen over storage duration. This enables better root cause analysis and quality decision-making across the product lifecycle.



## Background

### The Problem of Oxygen Ingress

Oxygen ingress after pasteurisation is a critical concern in beer packaging. Even trace levels of O<sub>2</sub> can:

- Degrade flavour and aroma compounds (oxidative staling)
- Reduce shelf life
- React with polyphenols, proteins, and other sensitive compounds
- Microbiological instability
- Loss of carbonation

After pasteurisation, bottles are sealed to preserve a CO<sub>2</sub>-rich environment (>99% CO<sub>2</sub> in headspace). However, oxygen ingress may occur due to

- Cap/crown sealing issues
- Leaks exacerbated by thermal expansion during pasteurisation
- Material permeability
- Handling and process variability

may allow atmospheric gas ingress, particularly oxygen and nitrogen.

Any oxygen ingress that occurs post-pasteurisation is critical because the product can no longer be thermally treated to mitigate its effects.

### Why Traditional O<sub>2</sub> Monitoring Falls Short

While direct measurement of headspace or dissolved O<sub>2</sub> can be effective shortly after packaging, it becomes increasingly misleading over time:

- Oxygen reacts chemically with beer components (e.g., aldehyde formation, polymerisation)
- These reactions consume the O<sub>2</sub>, reducing detectable levels
- This makes it difficult to determine whether ingress occurred weeks or months earlier

The **Orbisphere 6110** is a high-precision instrument designed for analysis of headspace gas composition in beverage packaging. Unlike traditional TPO or dissolved oxygen methods, headspace analysis:

- Directly assesses the gaseous environment above the liquid
- Identifies deviations in gas composition attributable to post-process ingress
- Allows process engineers and quality teams to pinpoint the timing and mechanism of contamination

### A New Approach: Residual Gas Fingerprinting

- Measures total gas volume and composition in the headspace
- Compares actual CO<sub>2</sub> purity against expected values
- Uses a formula to estimate non-CO<sub>2</sub> gas volume (i.e., foreign gas), which indirectly reflects historical oxygen ingress

This method can detect signs of **previous oxygen ingress even after the oxygen has reacted**, offering a powerful tool for process monitoring, QA diagnostics, and long-term storage evaluation.



## Method

### Equipment & Sample Handling

**Device:** Orbisphere 6110 Headspace Gas Analyzer

**Product Type:** Pasteurised beer in glass bottles

### Test Protocol

1. Analyse headspace gas composition immediately after pasteurisation.
2. Capture CO<sub>2</sub>%, O<sub>2</sub>%, and pressure data.
3. Apply foreign gas formula.
4. Evaluate test repeatability on identical batch samples.



Orbisphere 6110 Total O<sub>2</sub>/CO<sub>2</sub> Package Analyzer & Installation Kit

## Results and Analysis

### Results

Package Age	TPO (ppb)	CO <sub>2</sub> (True)	Vol. HS	Residual Gas (mL)	% Purity CO <sub>2</sub>
1 day	52.69	5.37	18.2	2.2	98.3
1 day	65.9	5.35	18.78	2.5	98.2
1 week	73.89	5.39	18.85	4	97.4
1 week	63.26	5.41	18.22	3.9	97.6
2 weeks	119.74	5.34	18.75	4.4	97.2
2 weeks	92.59	5.33	18.54	4.5	97.1
3 weeks	108.53	5.41	18.09	3.7	97.6
3 weeks	83.18	5.45	17.54	3.7	97.5
4 weeks	93.84	5.33	17.53	3.7	97.5
4 weeks	82.45	5.35	17.2	3.9	97.4
3 months	87.23	5.33	18.36	4.6	97.2
3 months	58.27	5.32	20.13	4.2	97.5
3 months	133.51	5.44	17.039	4.8	96.9
3 months	115.24	5.43	19.51	4.9	97.2
6 months	95.48	5.09	18.96	5.9	96.5
6 months	227.7	5.05	18.48	6.4	96.1
12 months	440.45	4.87	21.47	6	96.7
12 months	63.06	5.08	17.85	4.9	96.8

### Analysis

#### Key Observations by Package Age

Age	Avg TPO (ppb)	Avg. CO <sub>2</sub> Purity (%)	Avg. Residual Gas (mL)
1 day	59.30	98.25	2.35
1 week	68.58	97.50	3.95
2 weeks	106.17	97.15	4.45
3 weeks	95.86	97.55	3.70
4 weeks	88.15	97.45	3.80
3 months	98.56	97.20	4.63
6 months	161.59	96.30	6.15
12 months	251.76	96.75	5.45

## Trends & Interpretations

### TPO (Total Package Oxygen):

- Increases significantly with time. Starting from ~59 ppb at 1 day and reaching >250 ppb at 12 months.
- Notable jumps at 2 weeks, 6 months, and 12 months, suggesting progressive oxygen ingress over time.

### CO<sub>2</sub> Purity:

- Starts at ~98.2% and steadily declines to ~96.7% by 12 months.
- Indicates gradual ingress of foreign gas, likely air (O<sub>2</sub> and N<sub>2</sub>).
- Drop is modest (~1.5%), but consistent, supporting the ingress hypothesis.

### Residual Gas Volume:

- Increases from ~2.35 mL at 1 day to ~6.15 mL at 6 months  
Suggests either:
  - More headspace gas due to liquid CO<sub>2</sub> loss or
  - Introduction of foreign gas into the package
  - Correlates with increasing TPO and decreasing CO<sub>2</sub> purity.



## Correlation Insights

Age	Avg. TPO (ppb)	Avg. CO <sub>2</sub> Purity (%)
TPO vs. Age	+0.91	Strong positive - oxygen ingress over time
CO <sub>2</sub> Purity vs. Age	-0.86	Strong negative - loss of CO <sub>2</sub> purity with time
Residual Gas vs. Age	+0.87	Strong positive - headspace gas increases as ingress occurs

These strong correlations support the foreign gas ingress hypothesis, likely due to:

- Cap permeability or microleakage
- Thermal cycling stresses on crown/cap integrity
- Storage conditions over time

## Conclusions

- O<sub>2</sub> ingress increases significantly over time, confirming the value of early and repeated headspace analysis.
- CO<sub>2</sub> purity drops consistently, validating its use as a diagnostic parameter.
- Residual gas volume increases are indicative of either CO<sub>2</sub> loss or external gas ingress.
- The Orbisphere 6110 is well-suited for tracking these parameters in a non-destructive and repeatable manner.
- Apply SPC methods to headspace data across aging studies for predictive quality management.



## Recommendations

- Establish CO<sub>2</sub> purity thresholds to flag early signs of ingress (e.g., alert if <97.0% within 2 weeks).
- Use TPO + CO<sub>2</sub> purity + residual gas volume for multi-parameter fingerprinting.
- Investigate closure integrity and bottle variability if TPO deviates early (before 2 weeks).
- Apply SPC methods to headspace data across aging studies for predictive quality management.

## Discussion

The Orbisphere 6110 enables precise identification of package-level gas composition post-pasteurisation. This makes it a valuable tool for:

- Verifying seal integrity
- Diagnosing sources of oxygen ingress
- Detecting closure or filler issues
- Supporting corrective actions and preventative maintenance

Because the technique focuses on headspace changes, it provides a sharper diagnostic view than dissolved oxygen or TPO methods alone.

## Recommendations

- Introduce headspace fingerprinting into routine QA for post-pasteurised lines.
- Use foreign gas formula to benchmark new packaging components or closure types.
- Integrate 6110 data with SPC systems to monitor trends and preempt issues.
- Expand testing to different beverage types and packaging formats.

## Conclusion

The Orbisphere 6110 enables highly sensitive headspace gas analysis, offering a novel “O<sub>2</sub> fingerprinting” approach to detecting foreign gas ingress in pasteurised beverage packaging. This supports enhanced quality control, faster root cause investigations, and ultimately, more stable and consistent products.

