



Pipet Filler, Safety Bulb

Продукт №:

1465100

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Unique cone adapter allows universal pipet filling. Simply squeeze rubber bulb to charge for filling, place pipet into cone adapter and insert into sample. Release the bulb to fill. Remove bulb and cover pipet end with finger for dispensing. Easy and safe. Autoclavable.

Спецификации

Какво има в кутията:

Includes autoclavable elastic wristband.

Обем / опаковка, размер:

1 each

Какво има в кутията

Includes autoclavable elastic wristband.