



## Lead reagent set, 5-150 µg/L Pb

Продукт №: 2375000 Опасен

ВGN Цена: Свържете се с нас

Не е показана дата на доставка







Acid soluble Lead in a potable water sample is first concentrated on a Fast Column Extractor. The Lead is then eluted from the Extractor and determined colorimetrically with an indicator.

## Спецификации

Shipping Method: Lead Reagent set

Бр. тестове: 20

Измерване, диапазон:  $5 - 150 \, \mu g/L \, Pb$  Инструмент: All except DR 900

DR/800 Series

 Метод, име:
 LeadTrak

 Метод, номер:
 8317

 Параметър:
 Lead