



Reagent set, aluminum, 50 mL sample, ECR method

Продукт №:

2443400

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



For the determination of Aluminium, Eriochrome Cyanine R Method. For 25ml sample.

For determination of Aluminum by the Eriochrome Cyanine R (ECR) Method. Hach Method 8326. Sample Size: 50 mL Range: 0.002 - 0.220 mg/L. Reagent set contains ECR Reagent Powder Pillows, Hexamethylenetetramine Buffer Reagent Powder Pillows and ECR Masking Reagent Solution for approximately 100 tests.

Спецификации

Бр. тестове:	100
Измерване, диапазон:	0 - 0,220 mg/L Al
Инструмент:	DR/3000, DR/2000
Метод, име:	Eriochrome Cyanine R
Метод, номер:	8326