



**Be Right™**



## Apparatus set, cold vapor mercury

Продукт №: 2674400  
BGN Цена: Свържете се с нас  
Повикване за дата на доставка

Required accessory for the determination of Mercury with Mercury Reagent Set.

For determination of Mercury by the Cold Vapor Mercury Concentration Method. Hach Method 10065. Range: 0.1 - 2.5  $\mu\text{g/L}$  Hg. Method requires Cold Vapor Mercury reagents and other accessories.

---

### Спецификации

Метод, име: Cold Vapor Concentration