



Citrate Buffer Powder Pillows

Продукт №:

2107669

Опасен

BGN Цена:

Свържете се с нас

Налично



Powder reagent in sealed packages for the determination of Manganese, cold periodate method. For 10 mL sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required).
Hach method 8034. Sample size: 10 mL Pack of 100 Powder Pillows.

Спецификации

Бр. тестове:	100
Метод, име:	Periodate Oxidation
Обем / опаковка, размер:	100 /pk
Платформа:	Прахови вложки