



Sodium periodate powder pillows for manganese, 25 mL, pk/100

Продукт №:

98499

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Powder reagent in sealed package for the determination of Manganese, cold periodate method, and Iron titrations. For 25 ml sample.

For determination of manganese by the periodate oxidation method. EPA approved for reporting wastewater analyses (digestion required).
Hach method 8034. Sample size: 25 mL Pack of 100 Powder Pillows.

Спецификации

Бр. тестове:	100
Метод, име:	Periodate Oxidation
Обем / опаковка, размер:	100 /pk
Платформа:	Прахови вложки