



## Chlorophosphonazo Indicator Solution, 500 mL

Продукт №:

2589549

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



For the determination of Ultra Low Range Hardness, 4-1000  $\mu\text{g/l}$   $\text{CaCO}_3$ .

For ultra low range hardness determination by the Chlorophosphonazo method. Hach Method 8374. Range: 4 - 1000  $\mu\text{g/L}$  Ca & Mg as  $\text{CaCO}_3$ . Ideal for use with Hach's Pour-Thru Cell. Requires CDTA Reagent, Cat. No. 2589636. 500 mL bottle

### Convenience

Rapid Liquid methods utilize liquid reagents, bottletop dispensers and pour-thru cells to provide simplified testing.

### Performance

Since the pour-thru cell provides a stationary optical system, Rapid Liquid methods typically provide better precision than methods that require movement of sample cells.

### Efficiency

Rapid Liquid methods are ideal for high-volume testing.

---

## Спецификации

Бр. тестове:	250
Измерване, диапазон:	4 - 1000 $\mu\text{g/L}$ $\text{CaCO}_3$
Метод, име:	Chlorophosphonazo
Метод, номер:	8374
Обем / опаковка, размер:	ml