



Reagent powder pillows, ChromaVer[®] 3, chromium, , 5 or 10 mL, pk/1000

Продукт №:

1271028

Опасен

BGN Цена:

Свържете се с нас

Налично



For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Спецификации

EPA compliant:	Да
Бр. тестове:	1000
Измерване, диапазон:	0,010 - 0,700 mg/L Cr ⁶⁺
Инструмент:	All except DR/3000, DR/2000, DR 900
Метод, име:	1,5-Diphenylcarbohydrazide
Метод, номер:	8023
Параметър:	Chromium
Платформа:	Прахови вложки