




CDTA Digital titrator cartridge, 0.800 M

Продукт №:

1440301

Опасен

EUR Цена:

Свържете се с нас

Налично



Cyclohexanediamine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.800 M. Digits = mg/L CaCO_3 in 100 mL. Digits x 0.001 = mmole.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Спецификации

Измерване, диапазон:	100 - 4000 mg/L CaCO_3
Метод, име:	CDTA Titration
Платформа:	Цифров титратор

Задължителни аксесоари

- Защитни ръкавици, размер L, сини, нитрил, без пудра, 50 броя (Item SM743L)