



Reagent set, manganese, HR, periodate oxidation, 25 mL

Продукт №:

2243200

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



For determination of Manganese by the Periodate Oxidation Method, Hach Method 8034. EPA Approved for wastewater reporting (digestion required) - Federal Register, 44, (116) 34193 (June 14, 1979). Range: 0.1 - 20.0 mg/l Mn. Sample Size: 25 mL. Set Contains: Buffer Powder Pillows (Citrate type) and Sodium Periodate Powder Pillows for 100 tests.

Спецификации

EPA compliant:	Да
Бр. тестове:	100
Измерване, диапазон:	0,1 - 20,0 mg/L Mn
Инструмент:	DR/3000, DR/2000
Метод, име:	Periodate Oxidation
Метод, номер:	8034
Обем на проба:	25 mL
Параметър:	Manganese
Платформа:	Прахови вложки