



Amino acid F dilution solvent, for rapid liquid silica, 475 mL

Продукт №:

2353011

Опасен

BGN Цена:

Свържете се с нас

Налично



For determination of Ultra Low Range silica by the Heteropoly Blue Rapid Liquid method. Hach method 8282 Range: 3 - 1000 µg/L SiO₂. Used with Amino Acid F Reagent Powder (2651155). Rapid Liquid[™] methods offer an economical way to analyze multiple samples and to measure trace concentrations. Requires Pour-Thru Cell.

Спецификации

Бр. тестове:	250
Метод, име:	Heteropoly Blue
Платформа:	Разтвор