



Total Alkalinity Test Strips, 0-240 mg/L, 50 tests

| | |
|------------|-------------------|
| Продукт №: | 2744850 |
| BGN Цена: | Свържете се с нас |
| Налично | |

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to drinking water, food, and beverage processors; boiler operators; and aquatic biologists. With Hach's alkalinity test strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site, without any measuring, set-up, clean-up, or chemical handling. Our alkalinity test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Contains 50 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Спецификации

| | |
|-------------------------------|--|
| Бр. тестове: | 50 |
| Единица: | 50 tests |
| Измерване, диапазон: | 0 - 240 mg/L |
| Модел на тестов комплект: | Test Strips |
| Най-малки инкременти, стъпки: | Steps: 0, 40, 80, 120, 180, 240 mg/L (ppm) |
| Параметър: | Total Alkalinity - Phenolphthalein |
| Платформа: | Тестови ленти |