



## Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Продукт №: 2745050  
BGN Цена: Свържете се с нас  
Доставка до 2 седмици

The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Chlorine exists in water as "free" chlorine (hypochlorite) or "combined" chlorine (chloramines). These chlorine test strips streamline water testing for field technicians, enabling swift analysis of multiple samples within minutes and facilitating on-site evaluations without the need for extensive preparation or cleanup. Widely utilized in laboratories globally, our chlorine test strips serve as invaluable tools for pre-test screening, swiftly detecting potential interferences with lab testing procedures. Each package contains 50 tests.

### Спецификации

Бележка под линия:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$ ; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Единица:	50 tests
Измерване, диапазон:	0 - 10 mg/L $\text{Cl}_2$
Модел на тестов комплект:	TS
Най-малки инкременти, стъпки:	Steps: 0, 0.5, 1, 2, 4, 10
Параметър:	Free and Total Chlorine
Платформа:	Тестови ленти