



Bag, Whirl-Pak* 5-1/2X9, pk/10

Продукт №:	1437297
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Спецификации

Капацитет:	720 mL
Обем / опаковка, размер:	10