



## Digital Titrator with plastic case

Продукт №:	1690001
BGN Цена:	Свържете се с нас
Налично	

Portable and robust digital Titrator for fast and precise titrations. High resolution (125µL/digit) and high accuracy ( $\pm 1\%$ ). Easy handling. Applications for more than 40 parameters with interchangeable cartridges.

Incl. plastic case, manual & accessories.

The Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field. Interchangeable Cartridges allow to change titrant in seconds. No cleaning or rinsing of burets is necessary. The number used for calculation of sample concentration is recorded directly from digital counter. For many methods, the displayed value represents the concentration directly.

**Accurate results - when and where you need them**

**Fast and Flexible**

**Versatile**

**Economical**

---

### Спецификации

Гаранция:	24 месеца
Обем на доставката:	Digital Titrator only
С включени аксесоари:	не
Тегло:	0,132 kg
Точност:	$\pm 1\%$ (Uncertainty of readings is 1 digit. Most samples require more than 100 digits)