



Bottle, BOD, 300 mL, 6/pk

Продукт №: 62106
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

Спецификации

Body Material:	Glass
Package size UOM:	6pk
Капацитет:	300 mL
Материал:	Borosilicate Glass
Метод, номер:	8043
Обем / опаковка, размер:	6 /pk