



## Bottle, BOD, 300 mL, 24/pk

Продукт №:

62124

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

---

### Спецификации

|                          |                    |
|--------------------------|--------------------|
| Body Material:           | Glass              |
| Package size UOM:        | 24pk               |
| Капацитет:               | 300 mL             |
| Материал:                | Borosilicate Glass |
| Обем / опаковка, размер: | 24 /pk             |