



Стандартен разтвор на натриев хлорид, 100 $\mu\text{S}/\text{cm}$, 500 mL

Продукт №:

2971849

BGN Цена:

Свържете се с нас

Налично

Conductivity standard solution, 100 $\mu\text{S}/\text{cm} \pm 1 \mu\text{S}/\text{cm}$ at 25°C, 49 mg/L ± 0.25 mg/L as Sodium Chloride (NaCl), (1 $\mu\text{S}/\text{cm} \sim 0.5$ mg/L TDS), 500 mL bottle.

Спецификации

Package settings:	Bottle
Концентрация:	49 mg/L ± 0.25 mg/L as Sodium Chloride (NaCl)
Обем / опаковка, размер:	500 mL
Описание:	Conductivity calibration solution
Параметър:	Conductivity
Платформа:	Разтвор
Специални бележки:	Observed shelf life after opening, following recommended storage and use: 3 months
Условия на съхранение:	10 - 25 °C (protect from light)