



**Be Right™**



## Pump, Vacuum, Hand-operated

Продукт №:

**1428300**

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

This unit is ideal for field use and is a real convenience in the lab. Pumping rate is 15 cc (0.50 oz) per stroke. Entire unit is corrosion-resistant, sealed, and self-lubricating. Nozzle fits standard 6 mm (1/4") I.D. tubing. Supplied with 61 cm (2') of clear plastic tubing. Choose with or without vacuum gauge.

---

### Спецификации

Обем / опаковка, размер:

1 each