



Be Right™



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Продукт №:

1472459

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Спецификации

Body Material:	Plastic
Package size UOM:	12pk
Капацитет:	500 mL
Материал:	PP
Обем / опаковка, размер:	12 /pk