



Bottle, BOD, black, 300 mL

Продукт №:	2086047
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Спецификации

Body Material:	Glass
Package size UOM:	1pk
Капацитет:	300 mL
Материал:	Borosilicate Glass
Обем / опаковка, размер:	1 each