



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 500 mL

Продукт №:

2089679

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Спецификации

Единица:	6/pk
Капацитет:	500 mL
Материал:	Glass Screw Cap