



## Flask, Erlenmeyer, polymethylpentene capacity 125 mL, 6/pk

|                                |                   |
|--------------------------------|-------------------|
| Продукт №:                     | 2089873           |
| BGN Цена:                      | Свържете се с нас |
| Не е показана дата на доставка |                   |

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

---

### Спецификации

|            |                   |
|------------|-------------------|
| Вид:       | Erlenmeyer        |
| Единица:   | 6/pk              |
| Капацитет: | 125 mL            |
| Материал:  | Polymethylpentene |