



Crucible, Gooch, 40 mL Capacity

Продукт №:	2098641
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Спецификации

Капацитет:	40 mL
Обем / опаковка, размер:	1 /pk