



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Продукт №: 2324663
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Спецификации

Body Material:	Polycarbonate
Капацитет:	1000 mL
Обем / опаковка, размер:	4 /pk