



pH Paper, 0 - 14 pH Range, 100 tests

Продукт №:

2601300

BGN Цена:

Свържете се с нас

Налично

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples. Measures 0 - 14 pH.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Спецификации

Бр. тестове:	100
Единица:	100 tests
Измерване, диапазон:	0 - 14 pH units
Метод, номер:	10056
Най-малки инкременти, стъпки:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Параметър:	pH
Платформа:	Тестови ленти