



Cylinder, graduated, 1000 mL, w/ handles

Продукт №: 2612953
BGN Цена: Свържете се с нас
Не е показана дата на доставка

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered” double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Спецификации

Вид:	Една метрична стойност
Измерване метод:	To Deliver/To Contain
Капацитет:	1000 mL
Материал на тялото (измервателния уред):	Polypropylene
Обем / опаковка, размер:	1 each
Сертификати:	None