



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Продукт №:	2613049
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Спецификации

Stem:	Blue Cap/Stem
Еденица:	6 pk
Капацитет:	500 mL
Материал:	LDPE
Обем / опаковка, размер:	6 pk
Цвят:	Transparent