



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Продукт №:	2634001
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (**Please Note:** These sample vials are not autoclaveable)

Спецификации

Диаметър:	43mm
Дължина:	93mm
Еденица:	200 /pk
Капацитет:	fill line at 100 mL (+/-2.5 mL)
Обем / опаковка, размер:	200 /pk
Специална функция:	Sample Vial, with Thiosulfate