



Thermometer, Enviro-Safe, Non-Mercury, -10 - 260 °C

Продукт №:	2635701
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Спецификации

Деления на термометъра:	1 °C
Дължина:	405 mm
Обем / опаковка, размер:	1 each
Температурен диапазон:	-10 °C - 260 °C
Тип на флуида в термометъра:	Non-Mercury