



Cylinder, graduated, mixing, certified, 50 mL

Продукт №: 2636341
BGN Цена: Свържете се с нас
Не е показана дата на доставка

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Спецификации

Вид:	Една метрична стойност, смесено
Измерване метод:	To Contain
Капацитет:	50 mL
Материал на тялото (измервателния уред):	Glass
Обем / опаковка, размер:	1 each
Сертификати:	Class A