



**Be Right™**



## **Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Blue Cap/ Stem, 6/pk**

**Продукт №:**

**2928202**

**BGN Цена:**

**Свържете се с нас**

Не е показана дата на доставка

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

---

### **Спецификации**

|                          |               |
|--------------------------|---------------|
| Stem:                    | Blue Cap/Stem |
| Еденица:                 | 6 pk          |
| Капацитет:               | 500 mL        |
| Материал:                | LDPE          |
| Обем / опаковка, размер: | 6 pk          |
| Цвят:                    | Clear         |