



Xenosep consumables kit for EPA method 1664A testing

Продукт №: 2943800
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Order the Xenosep Consumables Kit for a quick and easy way to replace your most used components.

Consumables kit includes: 24 Xenosep SPE Filters, 47mm diameter, 24 Transfer Pipets, 15 mL capacity per pipet, 24 Inline Columns with Sodium Sulfate, 10g per column.

Спецификации

Обем / опаковка, размер: 1 each