



## pH paper, 1.0 - 11.0 pH range

|                                |                   |
|--------------------------------|-------------------|
| Продукт №:                     | 39133             |
| BGN Цена:                      | Свържете се с нас |
| Не е показана дата на доставка |                   |

Provides simple, fast, economical method of pH measurement. Sharp color change for each 0.5 pH unit. Convenient for sample pH adjustments. Recommended for sewage water, not finished water. Five 4.6-meter (15') rolls per package.

---

### Спецификации

|                          |                                     |
|--------------------------|-------------------------------------|
| Измерване, диапазон:     | 1.0 - 11.0 pH units                 |
| Обем / опаковка, размер: | 5 4.6-meter (15') rolls per package |