



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 125 mm, circular, 100/pk, qualitative

Продукт №:

50657

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Спецификации

Диаметър:	125 mm
Материал:	Cellulose
Обем / опаковка, размер:	100 /pk