



Sanitary Mounting Kit, 3700 Series and pH Sensors

Продукт №:	MH018S8SZ
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

316 stainless steel sanitary 2" tee and heavy duty clamp. For use with the sanitary style 3700-series inductive/electrodeless conductivity sensors.

Спецификации

Материал:	Gasket: EPDM
Процесна връзка:	2 "
Специални бележки:	Composition: stainless steel 2" tee piece heavy-duty clamp special cap EPDM compound gasket