



Buffer Solution, pH 7.00 (NIST), color-coded yellow, 4L

Продукт №:

2283556

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

pH Buffer standard solutions for laboratory and field analyses

Quality Buffer Solution, colour-coded yellow for easy identification. NIST traceable. pH is 7.00 ± 0.02 at 25°C.

For pH meter standardization. Color-coded for easy identification. pH is ± 0.02 at 25°C. 4 Liters

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Спецификации

NIST проследимост:	Да
Package settings:	Bottle
pH стойност (и):	7.00
Кодиран цвят?:	Да
Концентрация:	7.00 ± 0.02 pH
Обем / опаковка, размер:	4l
Описание:	Buffer Solution (color-coded Yellow) (NIST)
Параметър:	pH
Специални бележки:	Observed shelf-life after opening, following recommended storage and use: 6 months
Условия на съхранение:	10 - 25 °C (protect from light)