



Nitrate Ionic Strength Adjustor Powder Pillows - Intellical, pk/100

Продукт №:

2984799

BGN Цена:

Свържете се с нас

Налично

Ionic Strength Adjustor for use with Intellical ISENO3181 or sensION+ 9663 nitrate ion selective electrodes. The ISA is critical for determining low level nitrate concentrations. Powder form does not dilute samples. Pack of 100 Powder Pillows.

Спецификации

Метод, име:

Ion Selective Electrode

Платформа:

Прахови вложки