



Be Right™



Total Alkalinity Test Strips, 0-240 mg/L, 50 tests

Продукт №:

2744850

BGN Цена:

Свържете се с нас

Налично

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to drinking water, food, and beverage processors; boiler operators; and aquatic biologists. Contains 50 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Спецификации

Бр. тестове:	50
Единица:	50 tests
Измерване, диапазон:	0 - 240 mg/L
Модел на тестов комплект:	Test Strips
Най-малки инкременти, стъпки:	Steps: 0, 40, 80, 120, 180, 240 mg/L (ppm)
Параметър:	Total Alkalinity - Phenolphthalein
Платформа:	Тестови ленти