



**Be Right™**



## **pH Test strips, 4 - 9 pH units, 50 tests**

**Продукт №:** 2745650  
**BGN Цена:** Свържете се с нас  
**Налично**

Most natural waters range from pH 4 to pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range. Contains 50 tests.

**Easy to use, disposable, and inexpensive**

**Great for quick evaluations of pH levels in the lab or in the field**

---

### **Спецификации**

Единица:	50 tests
Измерване, диапазон:	4 - 9 pH units
Модел на тестов комплект:	Test Strips
Най-малки инкременти, стъпки:	Steps: 4, 5, 6, 7, 8, 9 pH
Параметър:	pH, mid range
Платформа:	Тестови ленти