



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Продукт №:

2890200

BGN Цена:

Свържете се с нас

Налично

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Спецификации

Бележка под линия:

*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.

Единица:

100 tests

Измерване, диапазон:

0 - 600 mg/L Cl_2

Най-малки инкременти, стъпки:

Steps: 0, 25, 50, 100, 200, 400, 600

Параметър:

Free Chlorine, High Range

Платформа:

Тестови ленти