



Be Right™



Desiccator cabinet, non-vacuum

Продукт №: 2618100
BGN Цена: Свържете се с нас
Доставка до 2 седмици

Stackable, polystyrene, transparent from any side. Door has foam rubber gasket, two roller latches for airtight fit. Jar of silica gel included. Inside dimensions: 21 x 18 x 13 cm (8.25 x 7 x 5"). Door opening: 19 x 11 cm (7.5 x 4.4"). Desiccator Grease available and sold separately.

Спецификации

Обем / опаковка, размер: 1 each