



Filter, glass fiber, diameter 55 mm, preweighed standard methods TSS procedure, 100/pk

Продукт №: 2546120
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Спецификации

Диаметър:	55 mm
Единица:	100/pk
Материал:	Glass Microfiber
Обем / опаковка, размер:	100 /pk
Пори, размер:	1.5 µm