



HSE-Pb, Lead Sensor Pack for HSA1000 Analyzer

Продукт №:

5040100

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Sensor Pack for Lead includes: 10 HSE-Pb scanning electrodes, 10 PrepTab PT-Pb sample preparation tablets, 10 sample tubes, and crushing rods. 10 tests.

Спецификации

ЕРА compliant:	Да
Бр. тестове:	10
Измерване, диапазон:	2 - 100 µg/L
Инструмент:	SA-1100
Метод, номер:	1001