



Sampling bag, sterile with dechlorinating agent, 100 mL, 100 pcs

Продукт №:

2075333

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10mg sodium thiosulfate (will dechlorinate up to 10 mg Cl_2 in 100 mL sample) four 100mL fill-lines.

Спецификации

Капацитет: 177 mL ml Fill Line

Обем / опаковка, размер: 100