



Incubator, Water Bath, 220 Vac

Продукт №: 2616302

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Temperature calibration is easy. Double-wall construction of stainless steel with polymer-coated tank, air-jacketed heating eliminates hot spot. Recessed heating element prevents burnout. Powerful centrifugal pump provides strong circulating

capabilities but low turbidity. Capacity: 22.1 liters. Exterior dimensions: 35.6 x 72.4 x 24.1 cm (14 x 28.5 x 9.5"). Interior: 29.2 x 49.5 x 15.2 cm (11.5 x 19.5 x 6"). Temperature range: 5 °C above ambient to 70 °C; 100 °C covered.

Temperature uniformity: ± 0.1 without cover at 37 °C.

Air-jacketed heating

Eliminates hot spots for more uniform heating

Easy to clean interior

Stainless steel interior

Maintains temperature and accuracy

Temperature uniformity of ± 0.1 °C at 37 °C.

Recessed heating element

Prevents burnout, which extends instrument life

Спецификации

Изисквания към захранването (напрежение):	220 V
Какво има в кутията:	Includes: Water bath, Gable cover, power cord and instructional manual.
Обем / опаковка, размер:	1 each

Какво има в кутията

Includes: Water bath, Gable cover, power cord and instructional manual.