



Paddle tester, total aerobic bacteria, yeast and mold, 10 pieces

Продукт №:

2610810

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Hach Paddle Testers allow fast, easy screening of water samples, solid surfaces, and nonviscous liquids where unsanitary conditions may exist. The paddle is a double-sided slide attached to the vial cap. Each side of the slide is used to perform a separate test. The slide has a molded-in grid for easy colony counting without opening the optically clear, leakproof vial. Paddle Tester, Total Aerobic Bacteria/Yeast and Mold (Tryptone Glucose Extract Agar/Sabouraud Dextrose Agar), pk/10.

Save time and simplify lab and field testing with Hach prepared P/A media

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Report Total Coliform values with EPA accepted Presence-Absence bottles and media

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Simplify routine testing with self-contained, easy-to-use Paddle Testers and BART products

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Спецификации

Media Type:

Total Aerobic Bacteria, Yeast and Mold

Обем / опаковка, размер:

10 /pk