



MEL/MPN Portable Laboratory

Продукт №: 2569800

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

MEL/MPN Lab

MEL/MPN Total Coliform and E. Coli Laboratory. Includes: portable Incubator, portable UV lamp, and consumables for 50 tests. **Media is not included (see Detailed Description).** For wastewater applications, order dilution water 14305-98.

Спецификации

Бр. тестове:	100
Чувствителност:	<1,1 MPN/100