




Calibration Cylinder Kit, Models 1720 E, D and C

Продукт №:

4415600

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Formazin Calibration Kit for 1720C and 1720D/E Turbidimeter. Includes calibration cylinder, Tensette Pipette, and 500 mL 4000 NTU Formazin.

Спецификации

Какво има в кутията:

Calibration cylinder, Tensette pipette 1.0-10.0 mL, and 500 mL 4000 NTU Formazin Standard, Instructions.

Какво има в кутията

Calibration cylinder, Tensette pipette 1.0-10.0 mL, and 500 mL 4000 NTU Formazin Standard, Instructions.